LipoNasal®

Heuschupfen

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DE

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Areas of application:

LipoNasal[®] Heuschnupfen serves to mitigate and concomitantly treat symptoms of hay fever (seasonal allergic rhinitis) as well as to humidify and care for the nasal mucosa. It is also suited for perennial allergic rhinitis (e.g. caused by animal allergens, mold and house dust mites).

LipoNasal[®] Heuschnupfen supports and promotes the natural healing process in irritation and inflammation of the nasal mucosa.

Directions for use:

Spray LipoNasal[®] Heuschnupfen 1 – 2 times into each nostril as needed. The treatment may be repeated several times a day. For more severe cases, LipoNasal[®] Heuschnupfen may be used more often.

Tips for correct use:

- Clean the inside of the nose before using the spray
 - Bend your head forward at an angle of approx. 45° («Looking at your own feet while standing»)
 - Hold the nasal spray upright. Place the head of the spray into one nostril and direct it slightly to the side of the nose (do not spray onto the nasal septum)
 - Inhale slowly through the nose and deliver a complete jet spray
 - Breathe out through the mouth
 - Repeat the same procedure on the other side of the nose.
- Please wipe off the head of the spray after each application and replace the protective cap after use.

In order to avoid possible infections, the spray bottle should only used by the same person.

Contraindications:

LipoNasal[®] Heuschnupfen should not be used in the presence of hypersensitivity to one of the ingredients.

Notes:

The bottle is filled generously, because for technical reasons it cannot be completely emptied. It is guaranteed that the indicated quantity can be used.

Information about the function of the membranes:

Inhaled air passes through the nose, the mouth and the throat (upper respiratory tract) as well as via the larynx, the trachea and the airways of the lungs (bronchial tract) into the two lobes of the lungs (lower respiratory tract). From here, inhaled air is transported via the fine branches (bronchioles) to the pulmonary alveoli (alveoli). It is here that the actual gas exchange takes place via the alveoli and the fine blood vessels (capillaries), i.e. oxygen is absorbed and carbon dioxide is released. The optimal conditions for gas exchange are prepared by the mucous membranes lining the nose. The mucous membranes moisten the inhaled air, clean and warm it, remove any detrimental influences and transport foreign bodies and secretions out of the lungs using tiny hair-like projections (cilia). The entire lining of our respiratory tract therefore inevitably comes into contact with bacteria, hazardous substances from the environment, as well as dust and pollen. The undisturbed function of the lower respiratory tract is also initially dependent upon the upper respiratory tract doing its task properly. Diseases of the upper respiratory tract very often continue as a rule in the lower respiratory tract.

The movement of the cilia in the mucous membranes of the respiratory tract constitute the key functions for cleaning and warding off infections. Inhaled particles are transported out of the nose by the rhythmic movements of the cilia. The cilia are the driving force of this "conveyor belt", although the membrane layer itself is the actual "conveyor". The normal mucous membrane film is divided into an internal, low viscosity saline layer, in which the cilia pulsate, and an outer, discontinuous gel layer (also called the transport layer), which is transported upwards by the movement of the cilia. Both layers are separated by a thin layer, the so-called surfactant. The liposomes contained in LipoNasal[®] Heuschnupfen consist of phospholipids, which represent about 75 % of the phospholipids of natural nasal surfactant in healthy persons. The surfactant is supported and stabilised by the addition of these phospholipids, which at the same time support the natural defence mechanism of the nasal membranes.

Shelf-life and storage:

Store in a cool and dry place. Keep out of reach of children. Can be used up to 6 months after first use. Do not use after expiry date. When using for the first time, you may have to pump 3 or 4 times to fill the tube.

Date of last revision: 2013-11